

Superior performance | Advanced features



8 electrodes



No coding



Can test neonates, pregnant women, anemia people etc.



Wide operating temperature



0.5 µL blood sample



5s testing time



Insufficient sample detection



Strip ejector



900-test memory



Meal markers



Hypo warning



5 test alarms

Specification

Feature	Specification
Measurement range	10 to 600 mg/dL (0.6-33.3 mmol/L)
Result calibration	Plasma-equivalent
Sample	Fresh capillary whole blood
Sample volume	About 0.5 µL
Test time	About 5 seconds
Power source	One (1) CR 2032 3.0V coin cell battery
Battery life	12 months or approximately 1,000 tests
Memory	900 records with date and time
Automatic shutoff	2 minutes after last action
Weight	Approximately 44.2g (with batteries)
Operating temperature	5-45 °C
Operating relative humidity	10-90% (non-condensing)
Hematocrit range	20-70%
Data port	Micro USB



Catalog | Ordering Information

Item	Catalog No.	Contents
VivaChek™ Ino Starter Kit	VGM01-011	Meter 10 Test strips 10 Lancets Lancing device Control solution N (2mL)
VivaChek™ Ino Starter Kit II	VGM01-012	Meter 10 Test strips 10 Lancets Lancing device
VivaChek™ Ino Meter Only Kit	VGM01-013	Meter 10 Lancets Lancing device Control solution N (2mL)
VivaChek™ Ino Meter Only Kit II	VGM01-014	Meter Lancing device
VivaChek™ Ino Meter Only Kit III	VGM01-015	Meter
VivaChek™ Ino Meter Only Kit IV	VGM01-016	Meter Lancing device 10 Lancets
VivaChek™ Ino Test Strips 50s (25/vial)	VGS01-011	50 Test strips (25/vial)
VivaChek™ Ino Test Strips 25s (25/vial)	VGS01-012	25 Test strips (25/vial)
VivaChek™ Ino Test Strips 10s (10/vial)	VGS01-013	10 Test strips (10/vial)
VivaChek™ Ino Test Strips 100s (25/vial)	VGS01-014	100 Test strips (25/vial)
VivaChek™ Ino Test Strips 50s (individual)	VGS01-015	50 Test strips (individual)
VivaChek™ Ino Test Strips 25s (individual)	VGS01-016	25 Test strips (individual)
VivaChek™ Ino Test Strips 100s (individual)	VGS01-017	100 Test strips (individual)
VivaChek™ Ino Test Strips 50s (50/vial)	VGS01-018	50 Test Strips (50/vial)
VivaChek™ Ino Test Strips 200s (50/vial)	VGS01-019	200 Test strips (50/vial)
VivaChek™ Ino Test Strips 100s (50/vial)	VGS01-01A	100 Test strips (50/vial)
VivaChek™ Ino Control Solution	VGC01-011	Control solution L (2mL) Control solution N (2mL) Control solution H (2mL)
VivaChek™ Lancing Device	VGD01-001	Lancing device
VivaChek™ Lancets 50s	VGL01-002	50 Lancets
VivaChek™ Lancets 100s	VGL01-003	100 Lancets
VivaChek™ Safety Lancets 50s	VGL02-001	50 Safety lancets
VivaChek™ Safety Lancets 100s	VGL02-002	100 Safety lancets

VivaChek
 VivaChek Biotech (Hangzhou) Co., Ltd.
 Level 2, Block 2, 146 East Chaofeng Rd.,
 Yuhang Economy Development Zone,
 Hangzhou, 311100, China
 www.vivachek.com

Number: 1280029401
 Effective Date: 2021-03-29
 ©2021 VivaChek Biotech (Hangzhou) Co., Ltd.

VivaChek™ Ino
 Blood Glucose Monitoring System



Simple & Reliable
 Fulfills all the requirements
 of EN ISO 15197:2015

Ecomed Medical (Pty) Ltd
 APPROVED
 Lizl Visage
 Training Facilitator





Cutting-edge technologies

The system detects the hematocrit signal of blood sample then reduces the interference within $\pm 10\%$

HCT Correction

The system detects the temperature signal then compensates the effect within $\pm 5\%$

Temperature Compensation

Accurate Results

Auto-coding

The system detects the coding info of each test strip automatically

Unique Formulation

We use unique formulation in the test strip which is very stable in extreme environments


Unique Algorithm

After detecting signals of various common sources of interferences, the meter runs unique algorithm and gives highly accurate result

Advanced 8 electrodes



The 8 electrodes on the VivaChek™ Ino test strip work with the VivaChek™ Ino meter to deliver highly accurate results. Detailed functions of those electrodes are:

- ✓ Code calibration (auto-coding)
- ✓ Eliminating hematocrit interference
- ✓ Eliminating temperature interference
- ✓ Checking humidity exposure
- ✓ Checking if the sample is sufficient
- ✓ Checking if the sample is blood or control solution
- ✓ Checking for possible damage of the test strip



Proven accuracy

A recent independent clinical study conducted by IDK Germany* shows that **98.8%** of all results from VivaChek™ Ino blood glucose monitoring system are within $\pm 15\%$ or 15 mg/dL (0.83 mmol/L) of the YSI plasma results, which is fully complying with the new **EN ISO 15197:2015** standard.



*IDK: Institute of Diabetes » Gerhardt Katsch « Karlsruhe (www.diabetes-karlsruhe.de).

Risk level (CEG zone)	Risk to diabetic patient
A	No effect on clinic action
B	Altered clinical action-little or no effect on clinical outcome
C	Altered clinical action-likely to effect clinical outcome
D	Altered clinical action-could have significant medical risk
E	Altered clinical action-could have dangerous consequences

System accuracy results for glucose concentration < 5.55 mmol/L (100 mg/dL)

Within ± 0.28 mmol/L (5 mg/dL) 78/186 (41.9%)

Within ± 0.56 mmol/L (10 mg/dL) 148/186 (79.6%)

Within ± 0.83 mmol/L (15 mg/dL) 185/186 (99.5%)

System accuracy results for glucose concentration ≥ 5.55 mmol/L (100 mg/dL)

Within $\pm 5\%$ 233/414 (56.3%)


Within $\pm 10\%$ 360/414 (87.0%)

Within $\pm 15\%$ 408/414 (98.6%)

System accuracy results for glucose concentration between 2.28 mmol/L (40.9 mg/dL) and 29.4 mmol/L (528.3 mg/dL)

593/600 (98.8%)

Hematocrit effect study



Over and under estimation of glucose due to the hematocrit effect is a common source of error. Anemia falsely elevates while polycythemia and dehydration falsely depresses blood glucose values. The magnitude of this effect may vary by 4-30% for every 10% change in hematocrit, depending on the systems.^{1,2,3} According to the WHO, anemia affects about 24.8% of the population in the world, about 47.4% of pre-school children and 41.8% of pregnant women worldwide.⁴

With the advanced Hematocrit Correction technology, the VivaChek™ Ino blood glucose monitoring system could be used for testing patients with hematocrit level within 20-70%, including neonates with very high hematocrit level, as well as anemia people and pregnant women with very low hematocrit level.

References:

1. Consensus Statement on Self Monitoring of Blood Glucose, Diabetes Care 1987; 1: 95-99.
2. The effect of Haematocrit on reagent strip test for glucose, Diabet Med 1991; 18 (2):172-175.
3. Factors Affecting Blood Glucose Monitoring: Source of Errors in Measurement H.Gnsberg Barry, MD, Ph.D.J Diabetes Sci Technol. Jul 2009; 3 (4): 903-913.
4. Worldwide prevalence of anemia 1993-2005, WHO Global Database on Anaemia Geneva, World Health Organization, 2008.

Ecomed Medical (Pty) Ltd
APPROVED

Lizl Visage
Training Facilitator

VivaChek™